**CERAVEST** is a carbon-free phosphate-bonded precision casting investment for crown and bridgework using all dental alloys by quick or progressive heating procedures.

<b>PHYSICAL DATA</b> (100% liquid concentration)	
Mixing ratio (powder/liquid)	: 100g /24ml
Flow	: 13 cm
Processing time (22-24°C)	: 5 – 7 min
Setting time (22-24°C)	: 9 – 12 min
Total expansion (linear)	: 3.00%
Compressive strength	: 4Mpa

**POWDER /LIQUID RATIO**

Ring size	Powder	Liquid / dist. water
x 1	1 x 75g	18 ml
x 3	1 x 150g	36 ml
x 6	2 x 150g	72 ml
x 9	3 x 150g	108 ml

**EXPANSION**

**Ceravest powder is mixed with Ceravest mixing liquid.** When this liquid is diluted with distilled water, the expansion of the Ceravest can be controlled so that compensation can be made for the contraction of the particular casting alloy used. The higher the concentration of mixing liquid, the greater the total expansion of the investing compound.

Type of alloy		Liquid / water ratio (%)	Ring size Liquid water volumes(ml)			
			x 1	x 3	x 6	x 9
<b>CROWN AND BRIDGES ALLOY</b>	High precious casting alloy > 70% gold	25/75	4 / 14	9 / 27	18 / 54	27 / 81
	Semi precious < 55% gold and Pd– base casting alloy	30/70	5 / 13	11 / 25	22 / 50	32 / 76
<b>CERAMIC ALLOY</b>	High/semi precious and Pd – base alloy	50/50	9 / 9	18 / 18	36 / 36	54 / 54
	Non– precious alloy NiCr	60/40	11 / 7	22 / 14	43 / 29	65 / 43
	Non– precious alloy CoCr	70/30	13 / 5	25 / 11	50 / 22	76 / 32

**STORAGE**

Store the powder and liquid at a normal room temperature (23°C). If the liquid is exposed to temperatures at or below 5°C, the liquid will freeze and cannot be used thereafter.

**PACKAGES**

Ceravest - powder 6 kg (40 x 150 g), liquid 1000 ml bottle

**DIRECTIONS FOR USE**

**PREPARATION** -Wetting agents are not necessary, but if used ensure that the wax is completely dry before investing.

**Metal ring method:** Use a moistened casting liner.

**MIXING** - Pre-mix the powder and liquid thoroughly by hand using a spatula. Mix for 60 seconds under a vacuum.

**INVESTING** - Invest within 6 minutes of the start of mixing (at 23°C). Invest under low vibration. When the ring is full, cease vibration and do not touch the investment until set. Leave to set for 20 minutes from the start of mixing. Scrape the top surface of the investment ring with a sharp knife.

**Ringless method:** After the initial setting (15 min. at 23°C), remove the plastic and allow to set fully.

**HEATING PROCEDURES**

**Quick heating method in a preheated furnace at the final temperature:**

Moulds (with a metal mould ring) of sizes x1 to x6 can be heated up rapidly.

Heat up the furnace according to the type of alloy to be used:

700 – 750°Cfor gold alloys

800 - 850°Cfor ceramic alloys

900°Cfor non precious alloys

After 20 minutes setting, place the investment directly into the furnace at the end-temperature.

**Heating method in a preheated furnace at 600°C:**

We recommended that moulds without ring are placed in a preheated furnace at 600°C, just 20 minutes after mixing. After that, heat directly to the final temperature (700°C – 900°C).

Hold at the end-temperature for:x 1 ...40 min x 3...50 min x 6...60 min x 9...90 min

Heating rate	x 1	x 3	x 6	x 9
Ambient to 250°C	5°C/min			
Holding time at 250°C	30*	40*	50*	60*
250°C to 570°C	7°C/min			
Holding time at 570°C	20*	30*	40*	50*
570°C to end temperature	9°C/min			
Holding time at end temp	30*	40*	50*	60*

**CASTING** - Cast in the usual manner: centrifugal, vacuum-pressure casting, etc... Cast as soon as possible after removing the ring from the burnout furnace. Cast according to the processing instructions from the alloy manufacturer.

**COOLING** - Place the ring upside down. Allow to cool to room temperature.

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