*Universal Composite Restorative*

*Universal Composite Füllung*

# **Dynamic** *Plus*<sup>®</sup>

**PREBOND**

**ETCHING GEL**

*Inhalt / Contents / Índice / Table Des Matieres /  
İçindekiler / محتويات / содержание*

*Universal Composite Restorative*

*Universal Composite Füllung*

# Dynamic *Plus*<sup>®</sup>

PRESIDENT DENTAL GmbH

Kesselbodenstrasse 5, 85391 Allershausen, Germany  
Tel: +49 8166 389 9 820 • Fax: +49 8166 389 9 821  
[www.presidentdental.com](http://www.presidentdental.com) • [info@presidentdental.com](mailto:info@presidentdental.com)

## • DYNAMIC PLUS

BENUTZERHANDBUCH .....	1	■
USER GUIDE .....	4	■
GUÍA DEL USUARIO .....	7	■
MODE D'EMPLOI .....	10	■
KULLANMA KILAVUZU .....	13	■
دليل المستخدم .....	15	■
GUIDA UTENTE .....	18	■
ГИД ПОЛЬЗОВАТЕЛЯ .....	21	■

## • PREBOND

BENUTZERHANDBUCH .....	24	■
USER GUIDE .....	25	■
GUÍA DEL USUARIO .....	26	■
MODE D'EMPLOI .....	27	■
KULLANMA KILAVUZU .....	28	■
دليل المستخدم .....	29	■
GUIDA UTENTE .....	30	■
ГИД ПОЛЬЗОВАТЕЛЯ .....	31	■

## • ETCHING GEL

BENUTZERHANDBUCH .....	32	■
USER GUIDE .....	33	■
GUÍA DEL USUARIO .....	34	■
KULLANMA KILAVUZU .....	35	■
دليل المستخدم .....	36	■
GUIDA UTENTE .....	37	■
ГИД ПОЛЬЗОВАТЕЛЯ .....	38	■

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### Dynamic PLUS

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## User Guide

Resin-based dental restorative material Description Dynamic is a visible-light activated, radiopaque, hybrid composite material. It is use in both anterior and posterior restorations. designed for

### Composition

- Barium aluminosilicate (Average Particle Size:  $\leq 1\mu\text{m}$ )
  - Fumed silica (Average Particle Size:  $0.04\mu\text{m}$ )
  - Bis-GMA
  - Triethyleneglycol dimethacrylate
- Others
- The percentage by weight of total inorganic iller is ca.80%

### Indications

Dynamic is indicated for class I, II, V of posterior teeth, class III, IV, V of anterior teeth and cervical cavities(wedge shaped defect) involving root surfaces.

### Contra-Indications

Dynamic should not be used with patients who develop hypersensitivity to it or have a history of hypersensitivity to methacrylate monomer.

### Recommended Procedure

#### 1. Shade selection

Clean teeth with pumice and water to remove surface stains or extraneous plaque. Teeth are not monochromatic. Consider shade and restoration depth. It is desirable to choose the shade after mocking up. Alternatively, Vita1 Lumin<sup>®</sup> Vacuum shade guide may be used. Using of rubber dam is recom mended for isolation.

#### 2. Cavity preparation

##### - Anterior restoration

Using the conventional acid etching, prepare cavity for all Class III, Class IV, and Class V restoration.

##### - Posterior restoration

Prepare the cavity. No residual amalgam or other base material should be left on the internal surfaces of preparation that would interfere with light transmission and the hardening of the restorative material.

#### 3. Pulp Protection

In deep cavities cover the dentin close to the pulp with a minimum amount of calcium hydroxide liner leaving the rest of cavity surface free for bonding. Glass lonomer or other Eugenol-free base materials may be used, if wished.

#### 4. Placement of matrix

Use a matrix system, preferably a transparent one, with proper wedging for proximal contracts. Pre-wedging is advocated to achieve slight separation and facilitate optimal proximal contact.

#### 5. Enamel and dentin treatment

Follow the manufacturer's instructions regarding etching, prim ing, adhesive application and curing. It is recommended to use Etching Gel and Prebond or Prebond SE

## 6. Dispensing the composite

### • Syringe

Dispense the necessary amount of restorative material from the syringe onto the mixing pad by turning the handle slowly in a clockwise manner. Immediately replace syringe cap. If not used immediately, the dispensed material should be protected from light. Place restorative material into the cavity using instrument contained.

### • Single dose capsule

Insert capsule into dispenser and rotate to gain the proper angle of entrance into cavity. Extrude restorative directly into cavity, using a slow and steady pressure.

## 7. Placement

Place and light cure restorative in increments in 2.5mm levels or less. For permitting extension of composite beyond cavity margins, overfill the cavity slightly. Avoid intense light in the working field. 1. Vita is a registered trademark of Vita zahnfabrik, Bad Sackingen, Germany.

## 8. Curing

Expose each area of restoration surface to a visible light source. Hold the light guide tip as close to the restorative materials as possible during light exposure. The recommended exposure time and maximum increment thickness for each shade is shown below.

## 9. Finishing

Immediately after curing contour restoration surfaces with fine finishing bur or stone. Surfaces with fine finishing bur or stone. Carefully adjust occlusion by removing material with a ine polishing bur or stone.

### Storage

1. Do not store at elevated temperature or intense light.
2. Do not store material in proximity to Eugenol containing product.

### Precautions

1. For dental use only
2. It should be used by a dentist and dental professional. Observe the procedure and the pre cautions.
3. Keep it away from a reach of non-dental user including children, senior citizens and others.
4. If skin sensitization occurs or allergy to acrylate occurs, stop using.
5. Exposure to uncured resin should be avoided. If contacts skin, wash immediately with soap and water.
6. If accidental contact with eyes or prolonged contact with oral tissues occurs, flush immediately with a large amount of water.
7. If refrigerated, the material must be allowed to reach room temperature before use.
8. Use of protective gloves and a no-touch technique is recommended.
9. Protective device for the respiration tract is highly recom mended before using.
10. In case of non-polymerization, check the light intensity of the light curing unit and then contact the nearest dealer(s) or manufacturer directly regarding the same quality matters.
11. Do not use after expiry date.

Shade	Hickness	Exposure Time	
		1.000 m W/cm <sup>2</sup>	400m W/cm <sup>2</sup>
A1, A2, A3, A3.5, B1, B2, B3, C2, C3, I, P, E	2.5mm	10 sec	20 sec
A4, D3, A2O, A3	2.5mm	15 sec	30 sec

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## User Guide

Prebond is an adhesive with outstanding adhesion to dentin and etched enamel. Its adhesion improves between the composite and tooth surface.

### Indication:

Prebond is used to improve the adhesion of composites on enamel and dentin. Further, it improves the adhesion between glass-ionomer underfill (lining) cement and composite fillings ("sandwich technique").

### Application:

The cavity should be prepared according to the rules for adhesive filling technique. Depending on the situation, apply calcium hydroxide paste or an underfill (lining) cement in deep cavities. The enamel is etched, rinsed and then dried in accordance with the etch bonding technique.

The glass-ionomer cement filling should not be etched, only roughened. Once the prepared cavity is dry, 1-2 drops of light-curing Prebond are dripped onto the mixing pad (immediately reclose the bottle) and a thin layer is applied to the tooth using a brush or foam pellet. Excess material is removed using a dry brush or a pellet. The material is then cured using halogen light (15 sec.).

### Side effects:

Undesirable side-effects are extremely rare if the product is prepared and used properly. Should you become aware of any undesirable side-effects, including in cases of doubt, please notify us, giving as precise a description of the attendant circumstances and symptoms as possible.

### Working instructions:

Ensure that the working instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the accompanying circumstances.

### Contra-indications:

Do not use Prebond in the area of the pulpa.

### Storage/shelf-life:

3 years, store in a cool, dry place. Protect against light!

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## User Guide

### Description

Etching Gel is phosphoric acid formulation thickened with natural polymeric materials. Etching Gel contains 37% phosphoric acid and is designed as a semi-gel type. Etching Gel can be obtained easy and available etching effect. Also Etching Gel leaves no residues on the surface of etched tooth contrary to silica thickened etchants.

### Composition

- Phosphoric acid
- Purified water
- Thickener
- Colorant

### Indications

Etching the enamel and dentine for adhesive restorations.

### Contra-indications

Do not apply on soft tissues including oral mucosa, eyes and dental pulp tissue.

### Recommended Procedure

1. Isolate tooth and prepare the cavity in a conservative manner.
2. Equip the disposable tip after removing the cap.
3. Etch enamel and dentin with etching agent for 15 seconds. Rinse gel thoroughly and remove excess water from the preparation with a gentle stream of air for 1-2 seconds.

4. Dry etched surfaces with clean air. Do not blow saliva onto etched surfaces. If etched surface do not appear chalky White, repeat etching. Salivary contamination must be avoided.
5. Repeat steps 3-4 if desired.

### Storage

1. Close the syringe cap after use.
2. Do not store at elevated temperature or intense light.

### Precautions

1. For dental use only.
2. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
3. Keep it away from reach of non-dental users including children, senior citizens and others.
4. Avoid contacting with eye and skin
5. If contacts skin, wash immediately with soap and water.
6. Using protective gloves is recommended.
7. If accidental contact with eyes or prolonged oral tissues occurs, flush immediately with a large amount of water and consult a physician.
8. If ingestion, give a large amount of water or milk, take an antacid, seek medical attention.
9. Do not use after expiry date.

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Universal Hybrid Composite Restorative


# Dynamic Plus<sup>®</sup>

## Ordering Information




PRD.100.01.100 Starter Kit	5x4gr syringe (A1, A2, A3, B1, B2)	1x5 ml. Prebond (Bonding Agent)	1x3 ml. syringe etching gel
PRD.100.01.100 Regular Kit	8x4gr syringe (A1, A2, A3, A3.5, B2, B3, A2Q, I)	1x5 ml. Prebond (Bonding Agent)	2x3 ml. syringe etching gel

PRD.100.01.102	A1	PRD.100.01.110	C2
PRD.100.01.103	A2	PRD.100.01.111	C3
PRD.100.01.104	A3	PRD.100.01.112	D3
PRD.100.01.105	A3,5	PRD.100.01.113	OA2
PRD.100.01.106	A4	PRD.100.01.114	OA3
PRD.100.01.107	B1	PRD.100.01.115	I
PRD.100.01.108	B2	PRD.100.01.116	P
PRD.100.01.109	B3	PRD.100.01.117	E


Apply acid conditioner  
followed by adhesive  
or self-etching adhesive.  
Light cure


\* Light output > 800 mW/cm<sup>2</sup>

## ETCHING GEL


## Application

Start at enamel

Extend to dentine


Remove gel and  
excess water


**Do not desiccate**

## PREBOND


## Application

Apply ample amount,  
leave undisturbed

Remove solvent with  
air syringe

Light Cure

Prebond is a one component dental universal adhesive for bonding dentin/enamel using adhesive restoration technique.