EN)

NOVA-TRAY®

Description

Novatray is a light-curing tray material in the form of prefabricated, baseplates to make rigid, custom impression trays and bite-blocks, as well as for other indications. It can be cured in all ultraviolet and Halogen light-curing units.

Composition

- Urethane dimethacrylate (UDMA)
- Glass
- Photo-initiators
- Pigment

Indications

1. Individual tray for use crown, bridge and total prosthetic technique.
2. Base plates for bite patterns and odontograms

Recommended Procedure

1. Mark margins on the gypsum model for individual tray. Use a space maintainer such as wax sheet if necessary
2. Fit the *Novo-Tray* on the gypsum model.
3. Cut the excess materials along a boundary line. Removed materials can be used for a making handle.
4. Put the gypsum model with tray resin into a UV light curing unit and polymerize for 3 to 5 minutes. Then remove the tray resin from the gypsum model, turn it over and expose the reverse side to the light. A longer light exposure can improve the physical properties of the material.
5. Trim the rough surfaces with stone or burs and bore a hole in tray as many as you need.


6. Remove any residues of wax and as light inhibition layer formed on the surface of the tray resin caused by the presence of oxygen of the air that can be removed with a tissue soaked in alcohol or acetone.

Storage

1. Store at dry and cool place.
2. Do not store at elevated temperature or intense light.
3. Put the lid back on after using it.

Precautions

1. For dental use only.
2. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
3. Keep it away from a reach of non-dental users including children, senior citizens and others.
4. Do not use *Novo-Tray* to a patient having inflammation or allergy against acrylate.
5. Exposure to uncured resin should be avoided. Put on the protective gloves using *Novo-Tray*.
6. If contacts skin, wash immediately with soap and water.
7. If accidental contact with eyes or prolonged contact with oral tissues occurs, flush immediately with a large amount of water.
8. Do not use after expiry date.

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