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Direction For Use

Cold Cure
Powder & Liquid**Instructions for use:****Cold Cure**

Mix Cold Cure Powder and Liquid to form a selfcure (autopolymerising) acrylic resin. Conforms to BS EN ISO 20795-1 Type II Class I For repairs, indirect relines and additions to removable prosthetic appliances. Will bond to acrylic teeth when cured by hydroflask.

Proportions and mixing

Dispense 10 drops of liquid into a mixing vessel and sprinkle powder onto the liquid whilst vibrating the dish until all the liquid is absorbed. Shake off any excess powder. Add 5 drops of liquid and stir slowly with a metal spatula for 10 seconds (45 seconds if room temperature is below 18°C) The mix should now be fluid, giving continuous flow when poured. Apply immediately.

Repairs

Cast the plaster model as usual. Chamfer the broken edges of the denture. If necessary protect both the fitting surface of the denture and the teeth adjacent to the repair with petroleum jelly. Coat model with plaster separating solution. Seat the denture on the model and secure with Sticky Wax if necessary. Flow the fluid mix into the fracture. (Do not use after mix ceases to flow freely)

Laboratory Relining

Cast the plaster model as usual. Produce locating matrix for denture in the usual way. Reduce fitting the surface of denture. Coat model with plaster separating solution. Make an appropriate size mix. Spread onto reline area of denture. When surface dulls locate on model using matrix. Seat with firm pressure for 10 seconds. Secure firmly and place in Hydroflask.

Other Applications

Replacing and adding teeth, addition of clasps.

Curing

For best results cure in hydroflask at 40°C for 10-20 minutes at 30psi Remove carefully, trim and polish. Will also self cure in 20 minutes.

Note: Prior to curing, ensure model has not been allowed to dry out as trapped air may affect the quality of acrylic.

Negative interactions with Dental Materials

None noted

Contraindications and side effects

Curing involves an exothermic reaction so care must be taken if it is handled whilst setting. Cured material may cause allergic reaction in those susceptible Methyl Methacrylate. May cause sensitisation by skin contact.

Precautions

For use by the dental professional only

Powder

Replace cap after use, Do not inhale dust, Store between 16-25°C

Liquit

Mixture based Methyl Methacrylate Monomer

Risk Phrases

Highly Flammable Flash point below 11°C Harmful by inhalation, contact with skin and if swallowed. Danger of cumulative effects. Irritating to eyes, respiratory system and skin. May cause sensitisation contact.

Safety Phrases

Keep container in well ventilated place. Keep away from sources of ignition-no smoking. Do not empty in to drains. Take precautionary measures against static discharge.

Avoid prolonged exposure to sunlight. in case of contact with eyes and skin: Flush with copious volumes of water. Seek medical attention. Store below 25°C. Replace cap immediately after use

Powder	Pro Cryla Powder
Liquid	Pro Cryla Liquid
Mixing Ratio	Powder: 21gr. / Liquid: 10 ml.
Dough Time	20 - 30 Minuten
Work Time	15 - 20 Minuten

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