DE

Zenit Color Composite®

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
EN

Zenit Color Composite®

Light-curing stains on composite basis Instructions for use

Range of application

Characterization of composite facings on crowns and bridges, individualization of tailor-made resin teeth as well as charac-terization of denture base resins in partial and full dentures. Available in ten different shades.

Instructions

The stains can be applied directly from the syringe. Alternatively, the stains can be put on a palette for mixing and applied with a thin brush. All shades are intermixable.

Zenit Color Composite Stains can also be mixed with composite pastes. Please avoid intro- ducing bubbles into the mixture and con-sider the instructions for the composite.

Apply **Zenit Color Composite** Stains in layers no thicker than 0.3 mm.

Polymerization

Zenit Color Composite Stains can be polymerized in all conventional light-curing systems. Before the final polymerization, the complete polymerization site may be covered with **Zenit Color** Cover. This avoids the formation of a new dispersion layer and facilitates the finalization of the object.

Use a light-curing system with an emission spectrum of 310 - 500 nm. Perfect physical properties can only be obtained with immaculate lamps. It is important to test the light intensity regularly, in line with the manufacturer's instructions.

Finishing
All silicone polishers (square-edge or knife-edge wheels, cylinder) as well as hard metal and diamond mills are suitable for finishing the material. The facing is polished with handpiece-mounted goat-hair brushes, woolen wheels and a polishing paste. Careful conditioning of the surface is essential if optimal results are to be achieved, thus virtually precluding deposits (nicotine, caffeine etc.) which could cause discoloration.


Guarantee
Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines. As our products are subject to continued development, we reserve the right to modify them.

Working range of the material

1- 3 minutes, depending on the light conditions.

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PRESIDENT DENTAL GmbH