Instructions For Use

General Description

Zenit Flex discs with RA Mandrels and strips are used by the Dental Surgeon for the finishing and polishing of Composite restorations.

The Zenit Flex Discs (Pop-On and Snap-On) impart an excellent polish to composites, microfilled materials, amalgams, and semiprecious and precious metals. The Zenit Flex Strip, because of the thin, flexible backing and center gapped abrasive can be successfully and safely used near the gingival, and in interproximal areas having tight contacts. Each Zenit Flex strips have two abrasive grits portion thus the finishing and polishing of a restoration requires fewer steps. Use of this all Zenit Flex systems will assure a final polish equal to or better than natural dentition.

Indications for Use

Zenit Flex Pop On and Snap On

Zenit Flex Finishing and Polishing Discs can be used to finish and polish anterior and posterior restorations. This system is especially useful to smooth the facial and lingual margins, substituting for or used in conjunction with rubber points. The inherent disc flexibility allows the disc to bend, mimicking the curvature of the tooth. These larger discs are useful when finishing and polishing indirect restorations in the dental laboratory.

Zenit Flex Strips

Zenit Flex Finishing and Polishing Strips because of the center gap abrasive can be used to finish and polish all interproximal restorative surfaces.

Precautions for Dental Personnel & Patients

1. Do not use any of the Zenit Flex Discs on dental handpieces capable of speeds greater than 35,000 rpms, or any high speed handpiece.
2. Use of a disc and mandrel at speeds greater than 35,000 rpms could cause the disc or mandrel to shatter which may result in injury.
3. Do not operate mandrel without disc attached.
4. Use of protective eyewear is recommended when operating handpieces.
5. As with any handpiece, to protect patients from injury, do not make contact with soft tissues.

The MSDS can be obtained from President Dental or contact your local subsidiary.

Instructions for Use

Before using the Zenit Flex Discs or Strips remove excess restorative material and contour the restoration to the desired shape using fine diamonds, or fluted carbide burs. The following procedure will produce a quality polish regardless of the Zenit Flex System used.

Zenit Flex Discs Instructions

Note:

- * Skipping a grit size in the polishing sequence may reduce the quality of the restoration's polish.
- * Zenit Flex products are color coded, from dark shades to light shades (finer grits).
- * Avoid touching composite with the mandrel or disc eyelet because a discoloration mark may occur. This mark can be removed by repetition of the polishing steps.
- * Between the usage of each grit, rinsing the restorative surface is recommended.
- * On average, time spent using each disc is approximately 15-20 seconds.

1. Place Zenit Flex disc on mandrel by gently pushing the eyelet portion onto the mandrel until the disc is secure. Abrasive side of disc can be facing up or down depending on need.
2. The polishing motion should be constant and one directional. A back and forth movement over the restorative enamel margin is not recommended
3. Use light pressure when polishing, let the discs do the work.
4. A Zenit Flex coarse grit disc can be used at about 10,000 rpm for gross reduction. Rinse and dry.
5. For final contouring and finishing utilize the Zenit Flex medium grit disc at about 10,000 rpms.
6. Next, change to Zenit Flex fine grit disc. The handpiece speed can be increased to about 30,000 rpms.
7. Finally, bring the restoration to a lustrous polish using Zenit Flex superfine grit disc at about 30,000 rpms for 15-20 seconds.
8. Wash away powder or debris from restorative surface.

Zenit Flex Strip Instructions

1. Insert gapped abrasive center of Zenit Flex finishing strip (blue/purple) between interproximal contact points.
2. Position the blue portion of the strip over the composite surface to be finished.
3. Firmly grasp both ends of the strip and draw the abrasive over the composite in a vigorous, seesaw motion.
4. Repeat procedure using the purple portion of the strip.
5. Discard used strip.
6. Repeat steps 1-5 with Zenit Flex polishing strips (pink/white) using the pink side first.

Storage and Use

Zenit Flex discs and strips are single use items. Properly dispose of used discs and strips in accordance with your healthcare facility's policy.

Mandrel sterilization

1. Sterilization Packaging: Mandrels may be loaded into dedicated bur trays and packaged. Standard dental sterilization packaging materials may be used such as cassettes, or pouches.
2. Autoclave: Listed below are four different steam sterilization cycle types that will effectively sterilize Zenit Flex disc mandrels: 15 minutes at 121°C/250°F wrapped
10 minutes at 121°C/250°F unwrapped 10 minutes at 132°C/270°F wrapped 3 minutes at 132°C/270°F unwrapped

The sterilization instructions above have been validated by PRESIDENT DENTAL as being capable of preparing these devices for multiple uses. It is the responsibility of the dental professional to ensure that desired results are achieved and the process is conducted using proper equipment, materials, and trained personnel. This requires validation and routine monitoring of the process. Any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

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