

INSTRUCTION FOR USE

Product Name: Zenit Brush - Pure Kolinsky Dental Composite Brush





Product Code:

- PRD.01.ZBR.0
- PRD.01.ZBR.000
- PRD.01.ZBR.1
- PRD.01.ZBR.SET.100

Product Description

Zenit Brush which is produced from pure Kolinsky enables you make perfect composite retouches.

Technical Specifications

	Product Code	Description	Ø	Length
1	PRD.01.ZBR.0	Glaze Brush #0 (Round)	1.2	7.5
				
2	PRD.01.ZBR.000	Glaze Brush #000 (Round)	1.2	5.5
				
3	PRD.01.ZBR.1	Opaque Brush #01 (Flat)	2.6	6.0
				
	PRD.01.ZBR.SET.100	Glaze Brush #0 (Round) + Glaze Brush #000 (Round) + Opaque Brush #01 (Flat)		
				

Instruction for Use

Indications for Dental Composite Brushes: Use the fine #1 brush for fine detailing and application of tints and opaques Use the #2 and #3 brushes to apply bonding and silanating agents and for cement removal prior to polymerization.

Package Details

- PRD.01.ZBR.0 : 1× Glaze Brush #0 (Round)
- PRD.01.ZBR.000 : 1× Glaze Brush #000 (Round)
- PRD.01.ZBR.1 : 1× Opaque Brush #01 (Flat)
- PRD.01.ZBR.SET.100 : 1× Glaze Brush #0 (Round)
1× Glaze Brush #000 (Round)
1× Opaque Brush #01 (Flat)

ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ

Наименование изделия: Zenit Brush — кисть из чистого колонкового ворса для стоматологических композитов

Код изделия:

- PRD.01.ZBR.0
- PRD.01.ZBR.000
- PRD.01.ZBR.1
- PRD.01.ZBR.SET.100

Описание изделия

Кисть Zenit Brush, изготовленная из чистого колонкового ворса, позволит сделать идеальное моделирование композитов.

Технические характеристики

	Код продукта	Описание	Ø	Длина
1	PRD.01.ZBR.0	Кисть для нанесения глазури №0 (круглая)	1,2	7,5
2	PRD.01.ZBR.000	Кисть для нанесения глазури №000 (круглая)	1,2	5,5
3	PRD.01.ZBR.1	Кисть для опака №01 (плоская)	2,6	6,0
	PRD.01.ZBR.SET.100	Кисть для нанесения глазури №0 (круглая) + кисть для нанесения глазури №000 (круглая) + кисть для опака №01 (плоская)		

Инструкция по применению

Показания к применению кистей для стоматологических композитов: используйте тонкую кисть №1 для тонкой детализации и нанесения оттенков и опаков

Используйте кисти №2 и №3 для нанесения бондинговых и силанизирующих агентов, а также для удаления цемента перед полимеризацией.

Информация об упаковке

- PRD.01.ZBR.0 : 1× Кисть для нанесения глазури №0 (круглая)
- PRD.01.ZBR.000 : 1× Кисть для нанесения глазури №000 (круглая)
- PRD.01.ZBR.1 : 1× Кисть для опака №01 (плоская)
- PRD.01.ZBR.SET.100 : 1× Кисть для нанесения глазури №0 (круглая)
1× Кисть для нанесения глазури №000 (круглая)
1× Кисть для опака №01 (плоская)